



Ringbio

African Swine Fever Virus ASFV Rapid Test Kits

African swine fever (ASF) is a devastating infectious disease of pigs, usually deadly. No vaccine exists to combat this virus. It does not affect humans nor does it affect other animal species other than pigs and wild boars. It can be transmitted either via direct animal contact or via dissemination of contaminated food (e.g. sausages or uncooked meat).

African swine fever is a OIE listed disease. This means it is a transmissible disease that has the potential for very serious and rapid spread, irrespective of national borders, is of serious socio-economic or public health consequence and is of major importance in the international trade of animals and animal products.

Here at Ringbio, the following products are available for your choices.

ASFV-Ab Rapid Test Kit

Product #: C30011

A rapid lateral flow immunoassay to detect ASFV antibody in 10min.

One drop of serum/plasma, with several drops of buffer

No instrument, no long time incubation

Package size: 50 tests.

ASFV-Ag Rapid Test Kit


Product #: C30021

A rapid lateral flow immunoassay to detect ASFV antibody in 10min.

One drop of serum/plasma, with several drops of buffer

No instrument, no long time incubation

Package size: 50 tests.



Ring Biotechnology Co., Ltd

E-mail: export@ringbio.com diego@nbgen.com Web: www.ringbio.com

Add: Building 3, Zhongtongtai TechnoPark, No. 11, Kechuang 14th St, Beijing 100176, CHINA

Tel: +86-10-56267496 Technical Support & Service: +86-18600362934