


Ringbio

Ring Biotechnology Co Ltd Company Profile & Product Information




The Company

RINGBIO is a pioneer biotechnology company in the food and animal diagnostic industry.

Located in the capital city of China, RINGBIO is a national high-tech company with more than 150 employees and laboratories more than 4500M²

RINGBIO is ISO9001 certified manufacturer. The ISO certificate scope covers
Production and Sale of Food Safety Rapid Test Kit/ELISA Kit,
Production and Development of Biological Rapid Test Kit/ELISA Kit for Scientific Research

NBGen,  NBGen are trademarks of RINGBIO.



Company Glance

RINGBIO has a full-set ability from design, research & development, production and sales.

From sequence to antigen, antibody to ELISA kits and lateral flow test strips, from DNA/RNA to PCR detection systems, we do all innovations by ourselves.



Dairy Testing

RINGBIO develops Rapid Test Kits and ELISA Kits for food safety.
In the field of dairy testing, RINGBIO is your reliable partner

Antibiotic Residue Tests

Mycotoxin Rapid Tests

Nutrients Rapid Tests

Milk Adulteration Tests

Other Tests



All dairy tests are now complied with **CE**



Dairy Testing

All dairy tests are 96T/kit
Combo test detects 2 drugs
TriTest detects 3 drugs
QuaTest detects 4 drugs




Antibiotic Residue Tests

Code	Product
100000	NBReader for Rapid Test Kit
100001	Beta Lactams Rapid Test
100002	Beta-Lactams + Tetracyclines, BT Combo Test
100005	Lincomycin Rapid Test
100006	Tetracyclines Rapid Test (High sensitivity)
100007	Tetracyclines Rapid Test (Low sensitivity)
100008	Sulfonamides Rapid Test
100009	Fluoroquinolones Rapid Test
100010	Neomycin Rapid Test
100011	Gentamicin Rapid Test
100012	Streptomycin Rapid Test
100013	Tylosin Rapid Test

Dairy Testing

All dairy tests are 96T/kit
 Combo test detects 2 drugs
 TriTest detects 3 drugs
 QuaTest detects 4 drugs



All dairy tests are now complied with 

Ultra Sensitive Beta/Tetra Combo Test Kit

Based on the next generation receptor immunochromatography to detect penicillins, cephalosporins and tetracyclines in one single test kit. All detection limits are ultra sensitive.

 **3/30 kinds**
detected drugs

 **37 °C**
incubation

 **8-10min**
detected time

 **EU MRL**
detected limit

 **Reader**
result visible

 **Cost effective**
amazing low cost

Drug	LOD(μg/kg)	MRL
Penicillin G	1-1.5	4
Ampicillin	2-3	4
Amoxicillin	2-3	4
Oxacillin	2-3	30
Cloxacillin	3	30
Dicloxacillin	2	30
Nafcillin	6-8	30
Procaine P	2	4

Drug	LOD(μg/kg)	MRL
Cefacetrile	8-10	125
Cefapirin	3-4	60
Cefazolin	18-20	50
Ceftiofur	40	100
Cefaquinome	5-6	20
Cefalonium	2	20
Cefaoperazone	2	50
Cefalexin	25-300	100

Drug	LOD(μg/kg)	MRL
Tetracycline	3-4	100
Oxytetracycline	2-4	100
Doxycycline	2-4	100
Chlortetracycline	3-4	100

Drug	LOD(μg/kg)	MRL
Tetracycline	60-80	100
Oxytetracycline	80-100	100
Doxycycline	40-60	100
Chlortetracycline	40-60	100

Dairy Testing

All dairy tests are 96T/kit
Combo test detects 2 drugs
TriTest detects 3 drugs
QuaTest detects 4 drugs



Antibiotic Residue Tests


Code	Product
100014	Chloramphenicol Rapid Test
100015	Beta/Tetra/Sulfa, BTS TriTest
100016	Beta/Tetra/AflaM1 TriTest BTM1 Rapid Test
100017	Beta/Tetra/Macrolides, BTM TriTest
100018	Beta/Tetra/CAP/Strep, BTCS QuaTest
100019	Beta/Tetra/Cefalexine, BTCex TriTest
100020	Linco/Quino/Tylo, LQT TriTest
100021	Beta-Lactams + Cefalexin, BCex Combo Test
100023	Macrolides Rapid Test
100024	Quino/Macro/Linco/Erythro, QMLE QuaTest
100025	Aminoglycosides Rapid Test, GNKS QuaTest

✂ RINGBIO has the ability to develop the test kits with different combinations you need. Contact us today to get your desired products.

Dairy Testing

All dairy tests are 96T/kit
Combo test detects 2 drugs
TriTest detects 3 drugs
QuaTest detects 4 drugs



All dairy tests are now complied with 

Mycotoxin Rapid Tests

Code	Product
------	---------

100004 Aflatoxin M1 Rapid Test

RINGBIO AflaM1 Rapid Test detects Aflatoxin M1 at several different levels, from 0.05ppb to 0.1, 0.2, 0.3, 0.4 and 0.5ppb, which meet EU, USA, Codex and many other MRLs.

Nutrients Rapid Tests

Code	Product
------	---------

E6001 Bovine Lactoferrin (bLF) ELISA Test

R6001 Biotin (Vitamin B7) Rapid Test Kit

R6002 Folic acid (Vitamin B9) Rapid Test Kit

R6003 Vitamin B12 Rapid Test Kit

Dairy Testing

All dairy tests are 96T/kit
Combo test detects 2 drugs
TriTest detects 3 drugs
QuaTest detects 4 drugs



Milk Adulteration Tests

400003 Cow Milk Adulteration Test

to detect cow's milk from goat/sheep milk, detection limit starts from 0.05%.

applicable for liquid milk, cheese.

particularly for dairy manufacturers to control milk adulteration.

Other Rapid Tests

Code	Product
100003	Beta-Lactamase Rapid Test
400002	Melamine Rapid Test
400001	Aflatoxin B1 Rapid Test
400009	Benzoid Acid Rapid Test

Dairy Testing

NBReader[®] Rapid Test Strip Reader

Product Code 100000

NBReader[®] is a new Point-of-Care analysing system. It is based on scanning the image of the test strip to produce digitalized data for all rapid tests.

It is used for rapid test cards or test strips based on immunochromatography tech. With this reader, users can obtain rapidness, high sensitivity, repeatable detection and avoiding mistakes by visual interpretation.

Technical details

Weight	3.0kg
Dimension	320mm×195mm×156mm
Wavelength	525±30nm
Precision	C.V. ≤3%
Display	6.2Inch LCD
Printer	Built-in thermal printer
Memory	4GB built-in / flash drive



Our Business

RINGBIO now connects to world dairy customers through a wide network of distributors and our own business team. For the last 6 years, we have established strong relationships with customers from more than 20 countries. The more connections we have, the more responsibility we carry.




To join the RINGBIO family, please contact us by email to export@ringbio.com. Always visit our website www.ringbio.com for news and information updates. Contact us at anytime, anywhere. Your attention is our priority.

How to Reach Us

RINGBIO is located in the Business Development Zone of Beijing, which is about 50Km from Beijing International Airport, and about 25Km from Beijing Railway Station.

Taking Taxi to RINGBIO from the airport will cost about 150CNY, which is about 25USD.

Schedule a visit and meet our team.
Welcome to RINGBIO



Ringbio



Ring Biotechnology Co Ltd

Add: Building 3, Zhongtongtai TechnoPark

Kechuang 14St, Beijing 100176, China

Tel: +86-10-56267496 / +86-10-56242030

Web: www.ringbio.com

Email: info@ringbio.com