


Beta-Lactams & Tetracyclines Combo Rapid Test Kit**BT Combo Test – MRL****Product Code: 102002****Product Information**

This product is designed to detect beta-lactams and tetracyclines antibiotic residues in raw milk. The detection limit meets European MRLs. The total time cost is 10min, no incubation, no special instrument is required. Very suitable for field testing or farm testing.

Technical Information

β-lactams	MRL(μg/L)	LOD(μg/L)	Cephalosporins	MRL(μg/L)	LOD(μg/L)
Penicillin G	4	2-4	Cefquinome	20	15-20
Ampicillin	4	3-4	Cefacetile	125	100
Amoxicillin	4	3-4	Cefalonium	20	18-20
Oxacillin	30	6-8	Cefazolin	50	50
Cloxacillin	30	6-8	Cefoperazone	50	40-50
Dicloxacillin	30	6-8	Cephapirin	60	50-60
Nafcillin	30	20-30	Ceftiofur	100	90-100
Tetracyclines	MRL(μg/L)	LOD(μg/L)	Tetracyclines	MRL(μg/L)	LOD(μg/L)
Tetracycline	100	40-60	Doxycycline	100	40-60
Oxytetracycline	100	80-100	Chlortetracycline	100	100

Schematic Steps**Package Information**

96 tests in one package.

In each package, there are 12 plastic tubes, in each tube, 8 test strips and 8 microwell reagents are included. A pack of 100pcs plastic pipette or micropipette tips also included.

Ordering Information**Ring Biotechnology Co Ltd**

Add: Building 3, Zhongtongtai TechnoPark, No 11, Kechuang 14St, Beijing 100176, China

Tel: +86-1056267496 / 56242030 Skype: ringbio

Web: www.ringbio.com Email: info@ringbio.com